


# System components in the CTL



# System components in the CTL

Some basic principles:

- Components **do not interact** with each other
- Each component manages **only its own** characterizing **properties**
- 4 principle types: **gantries**, **detectors**, **sources**, and beam modifiers
- Each type has a **general interface**
- Sub-classes are used to make **trade-offs** between flexibility and user convenience

*Base class  
(abstract)*

- Defines all required interfaces used throughout the CTL.
- Includes pure virtual methods that need to be implemented by derived classes.
- Can contain other virtual methods that can be overridden if required.

**non-instantiable**

*Generic*

- a generic, non-restrictive representation of the component type.
- Implements (pure) virtual methods from base class.
- Provides flexible but complex manipulation (setter).

**full flexibility  
complex interface**

*Specific*


- specific, “task-related” representations of the component type.
- Implement (pure) virtual methods from base class.
- Provide simple & convenient manipulation (setter).

**restrictive usage  
convenient interface**

## DETECTOR


## GANTRY

## SOURCE




## Background: Geometry - coordinate systems

### WORLD COORDINATE SYSTEM (WCS)



### CT COORDINATE SYSTEM (CTS)



## Background: Geometry – mat::Location

### mat::Location

- stores the position of a component in WCS
- also contains a rotation matrix describing the orientation
  - used for different purposes
    - [detector modules] transformation from a single module to whole detector
    - [gantry: source] transformation from CTS to WCS
    - [gantry: detector] transformation from WCS to CTS


```
struct Location
{
    Vector3x1 position = Vector3x1(0.0);
    Matrix3x3 rotation = mat::eye<3>();

    Location() = default;
    Location(const Vector3x1& pos, const Matrix3x3& rot);

    QVariant toVariant() const;
    void fromVariant(const QVariant& variant);
};
```

```
(
    typedef mat::Matrix<3, 3> Matrix3x3;
    typedef mat::Matrix<3, 1> Vector3x1;
)
```

# GANTRY



# GANTRY

## AbstractGantry

- manages the location of a source and detector component in WCS
- does not know what kind of source and detector is “mounted”

```
// abstract interface
protected:virtual mat::Location nominalDetectorLocation() const = 0;
protected:virtual mat::Location nominalSourceLocation() const = 0;
```

- provide the unmodified locations of the detector and source component  
→ not intended to be used as typical getter (--> protected); instead:

```
// getter methods
mat::Location sourceLocation() const;
mat::Location detectorLocation() const;
```


```
// convenience getter
Vector3x1 sourcePosition() const;
Matrix3x3 sourceRotation() const;
Vector3x1 detectorPosition() const;
Matrix3x3 detectorRotation() const;
```

- return the final location of the components, ie. including all potential **displacements** (later)



# GANTRY

## AbstractGantry



# GANTRY

## AbstractGantry

- **Displacements:**  
additional geometric distortions applied to the components in addition to their nominal geometry

```
// setter methods  
void setDetectorDisplacement(const mat::Location& displacement);  
void setGantryDisplacement(const mat::Location& displacement);  
void setSourceDisplacement(const mat::Location& displacement);
```

→ can be used for several purposes, eg.:

- Miscalibration
- Motion
- Quarter-detector-shift

### SOURCE geometry

$$t_{\text{src}}^{\text{total}} = R_{\text{gantry}}^{\text{displ}} \cdot t_{\text{src}}^{\text{nominal}} + R_{\text{src}}^{\text{total}} \cdot t_{\text{src}}^{\text{displ}} + t_{\text{gantry}}^{\text{displ}}$$

$$R_{\text{src}}^{\text{total}} = R_{\text{gantry}}^{\text{displ}} \cdot R_{\text{src}}^{\text{nominal}} \cdot R_{\text{src}}^{\text{displ}}$$

### DETECTOR geometry

$$t_{\text{det}}^{\text{total}} = R_{\text{gantry}}^{\text{displ}} \cdot t_{\text{det}}^{\text{nominal}} + \left(R_{\text{det}}^{\text{total}}\right)^T \cdot t_{\text{det}}^{\text{displ}} + t_{\text{gantry}}^{\text{displ}}$$

$$R_{\text{det}}^{\text{total}} = \left(R_{\text{det}}^{\text{displ}}\right)^T \cdot R_{\text{det}}^{\text{nominal}} \cdot \left(R_{\text{gantry}}^{\text{displ}}\right)^T$$

# GANTRY

Generic

## GenericGantry

- Allows free specification of the location of detector and source component

```
// setter methods  
void setDetectorLocation(const mat::Location& location);  
void setSourceLocation(const mat::Location& location);
```

Specific

## TubularGantry


- Source and detector component move on circular orbit with fixed diameter
- Orbit can be tilted around x-axis
- Configuration is fully specified by:
  - (gantry) rotation angle
  - (gantry) tilt angle
  - (table) pitch position

```
// setter methods  
void setPitchPosition(double position);  
void setRotationAngle(double angle);  
void setTiltAngle(double angle);
```


# GANTRY

## TubularGantry

```
void setRotationAngle(double angle);
```



```
void setPitchPosition(double position);  
void setTiltAngle(double angle);
```




# GANTRY

## CarmGantry

- Source and detector component move as a union with a mutable distance called “C-arm span”
- Configuration is fully specified by:
  - C-arm span
  - Location of the source (ie. position in WCS and rotation matrix to get there)

```
// setter methods  
void setLocation(const mat::Location& location);  
void setCarmSpan(double span);
```

# DETECTOR



# DETECTOR

## AbstractDetector

- Contains a set of identical flat panel detector modules  
→ all modules must have same number of pixels and same pixel size (and same skew)
- Each module has its own location, specifying its geometry with respect to the detector system as a whole (see next slide)

```
// abstract interface  
public:virtual QVector<ModuleLocation> moduleLocations() const = 0;
```

- Two (optional) response models
  - Saturation model → response in “intensity” domain (intensity, counts, or extinction)
  - Spectral response model → response in spectral domain

```
// setter methods  
void setSaturationModel(AbstractDataModel* model,  
                        SaturationModelType type);  
void setSpectralResponseModel(AbstractDataModel* model);
```

```
bool hasSaturationModel() const;  
bool hasSpectralResponseModel() const;
```

*considered within:*


DetectorSaturationExtension

SpectralEffectsExtension

# DETECTOR

## AbstractDetector

```
// abstract interface  
public:virtual QVector<ModuleLocation> moduleLocations() const = 0;
```





# DETECTOR


## GenericDetector

- Free specification of the location of individual detector modules


```
// setter methods  
void setModuleLocations(QVector<ModuleLocation> moduleLocations);
```

- Also allows for dynamic (ie. in between views) change of pixel sizes and allows setting of a non-zero skew factor


```
void setPixelSize(const QSizeF& size);  
void setSkewCoefficient(double skewCoefficient);
```



*original*



*changed size*



*non-zero skew*

# DETECTOR

## FlatPanelDetector

- Simplest implementation
  - All parameters fixed from moment of construction
    - Fixed pixel size
    - Always skew factor 0
    - Always a **single module** with location 0 (ie. no translation, no rotation)
- no setters

## CylindricalDetector

- Places individual detector modules on the surface of a cylinder  
→ represents typical curved CT detectors
- All parameters fixed from moment of construction
  - Fixed cylinder radius
  - Fixed number of modules
  - Fixed pixel size
  - Always skew factor 0

# DETECTOR

## CylindricalDetector

- Different ways of construction possible
  - Angulation + spacing
  - Radius + fan angle


```
static CylindricalDetector fromAngulationAndSpacing(const QSize& nbPixelPerModule,  
                                                    const QSizeF& pixelDimensions,  
                                                    uint nbDetectorModules,  
                                                    double angulationPerModule,  
                                                    double moduleSpacing,  
                                                    const QString& name = defaultName());
```

```
static CylindricalDetector fromRadiusAndFanAngle(const QSize& nbPixelPerModule,  
                                                  const QSizeF& pixelDimensions,  
                                                  uint nbDetectorModules,  
                                                  double radius,  
                                                  double fanAngle,  
                                                  const QString& name = defaultName());
```

# DETECTOR

## CylindricalDetector

- Different ways of construction possible
  - **Angulation + spacing**
  - Radius + fan angle




# DETECTOR

Specific

## CylindricalDetector

- Different ways of construction possible
  - Angulation + spacing
  - **Radius + fan angle**



# DETECTOR


## CylindricalDetector

- Different ways of construction possible
  - Angulation + spacing
  - Radius + fan angle
- several details on the configuration can be retrieved through various getters

```
// getter methods  
double angulationOfModule(uint module) const;  
double moduleSpacing() const;
```

```
// other methods  
double coneAngle() const;  
double curvatureRadius() const;  
double fanAngle() const;  
double rowCoverage() const;
```

## SOURCE



# SOURCE

## AbstractSource

- Source properties are characterized by
  - Energy range of emitted radiation (lowest to highest energy that is emitted)
  - Photon flux, ie. the total number of photons emitted to an area of  $1\text{cm}^2$  in one meter distance

```
// abstract interface
public:virtual EnergyRange nominalEnergyRange() const = 0;           [ typedef Range<float> EnergyRange; ]
protected:virtual double nominalPhotonFlux() const = 0;
```

- Similar to gantry classes, these are not intended to be used as (general purpose) getters, instead:

```
// getter methods
EnergyRange energyRange() const;
double photonFlux() const;
```

- Spectral properties are defined by a so-called spectrum model

```
// virtual methods
virtual IntervalDataSeries spectrum(uint nbSamples) const;
virtual uint spectrumDiscretizationHint() const;
virtual void setSpectrumModel(AbstractXraySpectrumModel* model);
```



# SOURCE

## AbstractSource

```
// virtual methods
virtual IntervalDataSeries spectrum(uint nbSamples) const;
virtual uint spectrumDiscretizationHint() const;
virtual void setSpectrumModel(AbstractXraySpectrumModel* model);
```

- setSpectrumModel() - setter for the spectrum model (can be restricted to certain types in reimplementations)
- spectrum() - returns a data series with values sampled from the spectrum model covering the full energy range with a given number of samples (note: relative photon counts!)
- spectrumDiscretizationHint() - returns a hint for a reasonable number of samples to use when sampling the spectrum (default 10)

- Some convenience getters:

```
float meanEnergy() const;
bool hasSpectrumModel() const; } self-explanatory
```

```
IntervalDataSeries spectrum(EnergyRange range, uint nbSamples) const;
```

- returns a data series with values sampled from the spectrum model covering only the specified energy range with a given number of samples

# SOURCE


## AbstractSource

- More general specifications:

```
// setter methods
void setFocalSpotSize(const QSizeF& size);
void setFocalSpotSize(double width, double height); // convenience alternative
void setFocalSpotPosition(const Vector3x1& position);
void setFocalSpotPosition(double x, double y, double z); // convenience alternative

void setFluxModifier(double modifier);
void setEnergyRangeRestriction(const EnergyRange& window);
```

- flux modifier: global (multiplicative) manipulation of the photon flux
- energy range restriction: restricts the energy range to a narrower interval
- focal spot specifications given in CT coordinate system



# SOURCE

## GenericSource

- Parameters can be set directly:

```
// setter methods
void setEnergyRange(const EnergyRange& range);
void setSpectrum(const IntervalDataSeries& spectrum, bool updateFlux = false);
void setPhotonFlux(double flux);
```

setSpectrum():

- data in spectrum is put into a TabulatedDataModel and used whenever a spectrum is requested from the component
- also updates the energy range to the smallest/largest value found in spectrum
- if updateFlux==true, integral over spectrum is used as new photon flux

→ useful, for example, to use externally generated spectra in the CTL

- another useful convenience feature:

```
static void setPhotonCountInSystem(SimpleCTSystem* system, double photonsPerPixel);
```

→ adjusts the photon flux in system such that each detector pixel in the current configuration receives (on average) photonsPerPixel photons (system must contain a GenericSource !)


# SOURCE

## XrayTube

- Represents a typical X-ray tube with parameters
  - Tube voltage → determines energy range and spectrum (rel. photon counts per energy)
  - Milliampere seconds (mAs) → determines photon flux

```
// setter methods  
void setTubeVoltage(double voltage);  
void setMilliampereSeconds(double mAs);
```

- Internally uses a spectrum model according to the TASMIP model [1]  
→ valid for tube voltages from 30 to 140 kV



# SOURCE

## — XrayLaser

- Source component that emits a photons of a single energy

```
// setter  
void setPhotonEnergy(double energy);  
void setRadiationOutput(double output);    → total energy emitted to 1cm2 in 1m distance (in mWs)
```

- Photon flux is computed from output power and photon energy

$\text{radiationOutput()} / (\text{ELEC\_VOLT} * \text{photonEnergy()} * 1.0\text{e}3)$

→ changing photon energy will change the flux, but keeps total emitted energy constant

- Energy range is  $[\text{photonEnergy}(), \text{photonEnergy}()]$
- `spectrumDiscretizationHint()` is 1